

MADE IN GERMANY FOR OVER 100 YEARS.
DEDICATED TO DIGNIFIED CARE.

VÖLKER

S 962-2 Hospital Bed





Home-like comfort even with increased care requirements.

The Völker S 962-2 Hospital Bed meets a broad range of requirements for use in hospitals. It allows patients to advance through the various steps of mobilisation, with support from caregivers — from sitting in the Comfort Chair position through to standing and getting out of the bed safely. A washable model, S 962-2 W, is also available.

By successfully combining form and function, Völker hospital beds can provide increased patient comfort while enabling caregivers to work safely and efficiently. The unique Healthcare Design from Völker ensures that all beds fulfil their medical purpose while creating a home-like atmosphere that encourages patient wellbeing and recovery.



Dignified care

Every person has a right to be treated with respect and dignity. Especially people in need of care. We think, develop and manufacture with this in mind.

Patented. The Völker siderail system.

The patented Völker siderail system is characterised by its exceptional versatility. It consists of separate elements that are split on both sides and can be adjusted to either one or two levels. This ensures that patients feel adequately protected from accidentally rolling out of the bed. As the siderails are highly versatile and robust, they can also assist mobilisation.

The side elements are directly integrated into the lying surface. The rails at the head and foot move with the bed as it is adjusted.



Comfortable. The Völker MiS® Micro-stimulation System.

The patented Völker MiS® Micro-stimulation System is based on neuroscience principles. It detects even minor movements made by the patient. The lying surface then feeds these back through the mattress and to the patient's body as tactile stimuli. This continuous feedback increases the patient's body perception, and promotes movement and blood circulation in the skin. Comfort is improved when the patient is sleeping or lying down, mobility is increased and the risk of bed sores is reduced.



Play it safe. With Völker Service.

Völker Service maximises the return on your investment in Völker beds. Our Service team will also take responsibility for a large amount of the work and requirements that arise from the German Medical Devices Act and other legislation.

All of our service technicians are qualified and

experienced testers who are very familiar with the medical devices they inspect. They will test and service your beds meticulously, observing all of the legal requirements.

Any defects are usually fixed immediately and therefore at optimum cost.

Völker Service additionally compiles detailed legal documentation and checks all records to ensure that they are complete.

We would be happy to discuss which Völker Service options would best suit your needs. Simply get in touch.

VÖLKER

MADE IN GERMANY FOR OVER 100 YEARS.
DEDICATED TO DIGNIFIED CARE.

Völker reserves the right to make changes to the design, specifications and models without prior notice. Products may vary slightly from the colours shown here. Subject to printing errors. All information is non-binding.

Not all products/options are available in all countries. For further information on our products or services, please contact your local Völker representative or visit our website.

©2018 Völker | All rights reserved

VÖLKER

www.voelker.de



**HEALTHCAREBEDS
AUSTRALIA**

info@healthcarebedsaust.com.au